

### Maximum fill additions to Infusion bags for doses prepared in clinical areas

The table below indicates the maximum fill additions which PASG recommends as guidance to be used by nursing staff working in clinical areas without the need to know the specific limits for an individual infusion bag brand. The information in the table below can safely be used for any infusion bag in use in the NHS listed below.

Bag size / mL	Maximum additional volume to be added to the infusion bag / mL
50	15
100	25
150	25
250	50
500	100
1000	100

The information in the table can be used for the following infusion container brands:-

Baxter Viaflo  
Baxter Viaflex  
Macopharma Easyflex N  
Macopharma Macoflex  
Macopharma Macoflex N  
Fresenius Freeflex  
B Braun Ecoflac

N.B. Pharmacy aseptic units have access to detailed information for these containers, and therefore may apply higher fill volumes when compounding aseptically prepared products.

Acknowledgment to Emer Sheridan for her contribution to this work.